


# StableVac™

## Fluid Management System

StableVac makes Fluid Management less invasive and provides superior distention, visualization and protection against extravasation during Arthroscopic procedures.



Patent Pending



### Advantages:

- Stable operative site for knee, shoulder and other joint procedures using gravity-fed irrigation infusion
- Superior Distention
- Superior Visualization
- Reduced risk of fluid extravasation due to elimination of infusion surge
- Reduced personnel time managing bags of irrigating fluids
- Less fluid leakage onto the floor
- **Cost Savings:**
  - Irrigation infusion pump and costly disposable tubing is eliminated
  - Fewer bags of saline used during procedure



### StableVac Ordering Information

Catalog No.	Description	Surgical Use	Packaged
<b>A4000-8</b>	StableVac Suction Tubing – 10 foot length	Arthroscopic procedures	10/box, Sterile
<b>A4206</b>	Dual Irrigation Set - 9 foot length	Dual Irrigation Infusion Set	10/box, Sterile
<b>A4500</b>	StableVac Combination Pack Includes StableVac and Dual Irrigation Set	Arthroscopic procedures	10/box, Sterile


\* Flow rate estimates are for reference purposes only



37 Shield, Irvine, CA 92618 USA  
 800.753.7400 toll free in U.S.  
 714.832.6300 Fax 714.832.0300  
 info@stablevac.net  
 www.stablevac.net



# StableVac Setup



## StableVac can be used in the following Arthroscopic procedures

### Large Space Arthroscopy



Shoulder Arthroscopy



Knee Arthroscopy

### Small Space Arthroscopy



Ankle Arthroscopy



Wrist Arthroscopy



Elbow Arthroscopy



Hip Arthroscopy